



November 17, 2016

## **NEW PRODUCTS**

## Native and Mass-Labelled PFAS Solution/Mixtures PFAC-24PAR & MPFAC-24ES

In response to customer requests for comprehensive mixtures of perfluorinated compounds, Wellington has prepared native and mass-labelled solution/mixtures containing commonly analyzed per- and poly-fluoroalkyl substances (PFASs). The native solution/mixture (PFAC-24PAR) contains 24 components (C<sub>4</sub>-C<sub>14</sub> PFCAs, C<sub>4</sub>-C<sub>10</sub> PFSAs, FOSA, N-MeFOSAA, N-EtFOSAA, as well as 4:2, 6:2, and 8:2FTS). To facilitate isomer specific quantification of PFHxS and PFOS, this solution was designed to include characterized mixtures of the linear and branched isomers of the C<sub>6</sub> and C<sub>8</sub> perfluoroalkanesulfonates. The mass-labelled solution/mixture (MPFAC-24ES) is a 19 component mixture containing suitable surrogates for the PFAC-24PAR.

Catalogue Number	Product (methanol)	Qty	Conc
PFAC-24PAR	Native PFAS Precision and Recovery	1.2 ml	2.0 μg/ml ea
	Standard Solution (24 components)		
MPFAC-24ES	Mass-Labelled PFAS Extraction	1.2 ml	1.0 μg/ml ea
	Standard Solution (19 components)		

These new solution/mixtures complement our existing line of mixed PFAS reference standards.

Catalogue Number	Product (methanol)	Qty	Conc
PFC-MXA	Native PFCA Solution/Mixture (C <sub>4</sub> -C <sub>14</sub> )	1.2 ml	2.0 μg/ml ea
PFS-MXA	Native PFSA Solution/Mixture (C <sub>4</sub> ,C <sub>6</sub> -C <sub>8</sub> ,C <sub>10</sub> )	1.2 ml	2.0 μg/ml ea
PFAC-MXA	Native PFCA/PFSA Solution/Mixture (10)	1.2 ml	5.0 μg/ml ea
PFAC-MXB	Native PFCA/PFSA Solution/Mixture (17)	1.2 ml	2.0 μg/ml ea
PFAC-MXC	Native PFCA/PFSA Solution/Mixture (21)	1.2 ml	2.0 μg/ml ea
MPFAC-MXA	Mass-Labelled PFCA/PFSA Solution/Mixture (9)	1.2 ml	2.0 μg/ml ea
MPFAC-C-ES	Mass-Labelled PFCA/PFSA Extraction Standard (13)	1.2 ml	2.0 μg/ml ea
MPFAC-C-IS	Mass-Labelled PFCA/PFSA Injection Standard (4)	1.2 ml	2.0 μg/ml ea
EPA-537IS	U.S. EPA Method 537 Internal Standard PDS (3)	1.2 ml	variable
EPA-537SS	U.S. EPA Method 537 Surrogate PDS (3)	1.2 ml	variable



Please contact your local distributor or info@well-labs.com for pricing and delivery.

Visit our website (www.well-labs.com) for a complete listing of our new products.