





NEW PRODUCTS

DECHLORANE PLUS®

Dechlorane Plus® (1,2,3,4,7,8,9,10,13,13,14,14-dodecachloro-1,4,4a,5,6,6a,7,10,10a,11,12,12a-dodecahydro-1,4,7,10-dimethanodibenzo[a,e]cyclooctene) is a chlorinated flame retardant which has been used commercially for over 40 years in the coating of electrical wire and cables, computer connectors and plastic roofing materials. The technical Dechlorane Plus® product is a mixture of two stereoisomers, the *syn*- and *anti*-conformers.



syn-Dechlorane Plus®



anti-Dechlorane Plus®

Despite its long commercial history, Dechlorane Plus® has only recently been detected in the environment. The potential for environmental translocation, bioaccumulation, persistency and toxicity of Dechlorane Plus® needs to be assessed in greater detail to understand its impact on humans and wildlife.

DECHLORANE PLUS®

Catalogue Number	Product (toluene solution)	Qty/Conc
s-DP	<i>syn</i> -Dechlorane Plus®	1.2 ml 50 µg/ml
a-DP	<i>anti</i> -Dechlorane Plus®	1.2 ml 50 µg/ml

To help further determine the levels of Dechlorane Plus® in the environment and the potential implications, **WELLINGTON** has now added reference standards for the two Dechlorane Plus® stereoisomers to its product line.

Please contact your local distributor or info@well-labs.com for pricing and delivery.

AND...Please continue to visit our website (www.well-labs.com) for a complete listing of new product releases.